



Shielded Fume Hood SFH 1200

- Designed to house 1200 mm wide fume hood type IF/ZK 120
- Outside dimensions are 1585 mm x 971 mm x 2600 mm (wxdxh)
- Internal dimensions of the fume hood 1174 mm x 630 mm x 700 (wxdxh)
- Sliding body shield 400 mm wide with lead glass 300 x 250 mm
- Ergonomically sloped front and body shield
- Worktop out of 316 L stainless steel
- Inside hood completely out of phenolic resin, and resistant against virtually every corrosive liquid, gas or vapour
- Lead shielding of 50 mm on table
- Lead shielding of 50 mm in side walls and back wall up to 1000 mm above table level
- Air outlet: 200 mm diameter
- Two Lighting units 15W exchangeable from the inside
- Electrical dual outlet included on the left and right outside of the hood
- Safety glass panel, vertically adjustable
- Outside smooth painted steel finish, easy to decontaminate
- Steel foot pads 10 mm thickness for weight distribution
- Power connection 230 V / 50 Hz
- Total weight approx. 2.800 kg