

MicroLite® Series

MicroLiteMX™ Radiation Protective Eyewear

AttenuTech®

Versatile Frame Style with Widest Rx Range

ScaNrad

PROTECTION AS

(+47)22 22 24 80 post@scanrad.no



(l-r) MicroLiteMX
Navy Blue (FB1000BL),
Matte Black (FB1000B),
Silver (FB1000S), and
Crimson (FB1000C)



- Unique “Spring Relief” temple bars provide secure, comfortable grip for a wide range of face and head characteristics
- Accommodates extraordinary range of positive and negative prescriptions, including progressive bifocal
- Super soft, ultra comfort nose piece
- Integrated Optibloc™ leaded side shields provide lateral protection against scatter radiation
- ProLiteHT™ high light transmission leaded lenses for excellent optical acuity
- Accessory bonus: protective case, microfiber bag, anti fog lens cleaning solution, adjustable ClickFit™ retaining cord.

Additional MicroLiteMX Features

- Attenuation: 95% radiation reduction at 125 KVP Direct Beam
- Pb Equivalent Lens: 0.75 mm (+/- .05) Pb glass
- Weight: 58 grams
- Available in plano, single vision Rx, bifocal Rx and progressive bifocal Rx
- Widest range of prescriptions, with fastest Rx turnaround in the industry
- Frame colors: navy blue, matte black, silver, or crimson
- One year warranty from date of sale



AttenuTech is the only company that designs and manufactures eyewear specifically for your protection in fluoroscopic environments. All AttenuTech® MicroLite® Series radiation protective eyewear is prescriptionable by design, with the fastest Rx turnaround in the industry.

Exclusive!
Only AttenuTech®
MicroLite® Series
Eyewear
has ProLiteHT™
Lenses

- Up to 20% more light transmission than all other leaded eyewear brands.
- Provides ultimate contrast and visual detail in the operating room.
- The lightest weight lenses available anywhere.